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





MD

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**EN****DRAINAGE SYSTEM  
TO BE USED  
FOR SUCTION ONLY****INTENDED USE**

A continuous suction drainage system of the surgical wound to prevent post-operative haematomas.

**DESCRIPTION**

MINI EMODREN is most suited to surgery of small parts of the hand or foot, micro-surgery and reconstructive cosmetic surgery for clinical purposes.

MINI EMODREN 50 ml is supplied with a 7 cm suction tube with a clamp, a manually detachable multi-diameter connection, and a drainage tube TD 14/22 (diameter 2.4 mm Ch 7) that adapts perfectly to every anatomical profile, including fingers. On request, the set can be provided with a 3.2 mm (Ch 10) diameter drainage tube.

MINI EMODREN 100 ml is available with drainage tubes TD 14/22 (diam mm 2.4 Ch 7), TD 14/32 (diam mm 3.2 Ch 10), TD 12/42 (diam mm 4.2 Ch 13) and TD 12/47 (diam mm 4.7 Ch 14).

Drainage tubes are produced without intentionally added BPA.

The choice of model on the basis of patient characteristics, surgery and the correct size is up to the physician. Device intended for adult patients.

The **P version**, for both the 50 ml and 100 ml models, is packaged with a disposable, atraumatic stainless steel guide needle already assembled on the drainage tube.

**INSTRUCTIONS FOR USE:**

- A. BELLWS**
- B. SUCTION TUBE CONNECTION**
- C. CLAMP**
- D. MULTI-DIAMETER CONNECTION**

1. Introduce the drainage tube into the surgical wound from the inside to the outside by means of its guide needle of the same diameter.

**Always employ strict asepsis.**

2. Fasten the drainage tube far from the wound with a single suture or adhesive tape. Connect the drainage tube projecting from the skin to the suction tube by means of the special connector, having first chosen the section of the right diameter. The various sections of the connection can be removed by gripping firmly and unscrewing (Fig. 4).

**3. Take care to shut the clamp.**

4. Compress the bellows and connect up the suction tube connection by screwing it into the hole on the top of the bellows.

**The clamp must be kept closed at this stage.**

5. Start suction by opening the clamp

**CONTRAINDICATIONS**

There are no known contraindications.

**WARNINGS**

- Re-use of single devices may involve the following risks:
  - deterioration of the performance of the devices;
  - infection or contamination by previous use on another patient.
- The device must only be used by qualified medical staff.
- After use, the device must be treated and disposed as a hospital waste, taking all the necessary precautions and as required by current provisions of law. Sharp parts must be disposed of in special rigid containers.
- The bellows is graduated to check the quantity of blood suctioned. The quantity in ml must be considered approximate.
- The use of original components is recommended. Medicalplastic declines all responsibility for harm caused to patients by the use of 3rd party products or use other than those specifically prescribed.
- The device is supplied sterile. Do not use after the expiry date or if the packaging has been opened or damaged.
- In case of a serious incident, please promptly notify the manufacturer and your Competent Authority.
- The device has not been specifically tested on pregnant and lactating women and paediatric patients. Any use on these populations should be evaluated by the physician.