



emodren® 500-700 ml

MD

	Ø mm 2-2,6 (Ch 6-8)	TD 14/22
	Ø mm 3-4,3 (Ch 9-13)	TD 14/32 - TD 12/42
	Ø mm 4,6-5,7-6,0 (Ch 14-17-18)	TD 12/47 - TD 12/57 - TD 12/60

Medicalplastic s.r.l.
Via Saverio Mercadante, 15 - 20124 Milano - Italia
tel. 026704424 - Email: info@medicalplastic.it - www.medicalplastic.it
MADE IN ITALY

500800

Rev.: 0 25/11/2024

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EN DRAINAGE SYSTEM TO BE USED FOR SUCTION ONLY

INTENDED USE

A continuous suction drainage system of the surgical wound to prevent post-operative haematomas.

DESCRIPTION

It is available with a 100 cm long suction tube, a manually detachable multi-diameter connection, a drainage tube, and a double Y-connection for an eventual second drainage tube.

The full range of Emodren products makes it possible to assemble different sets combining 500/700 ml bellows and 50 cm long drainage tubes of different diameters with perforations along 14 cm (single rows of diameters measuring 2.4 mm (Ch 7) and 3.2 mm (Ch 10) – double offset rows of diameters measuring 4.2 mm (Ch 13) – 4.7 mm (Ch 14) – 5.7 mm (Ch 17) – and 6.0 mm (Ch 18).

Drainage tubes are produced without intentionally added BPA.

The choice of model on the basis of patient characteristics, surgery and the correct size is up to the physician. Device intended for adult patients.

The **Base version** is without the drainage tube but can be supplied with two connections for silicone tubes.

Model P is packaged with a disposable Inox guide-needle already assembled on the drainage tube.

INSTRUCTIONS FOR USE:

- A. BELLWS
- B. SPRING
- C. SPRING FASTENING, UPPER
- D. SPRING FASTENING, LOWER
- E. SPRING EYELET
- F. SUCTION TUBE CONNECTION
- G. CLAMP
- H. SUSPENDER CLIP
- I. PLUG

1. Introduce the drainage tube into the surgical wound from the inside to the outside by means of its guide needle of the same diameter.

Always employ strict asepsis.

2. Fasten the drainage tube far from the wound with a single suture or adhesive tape.

Connect the drainage tube projecting from the skin to the suction tube by means of the special connector, having first chosen the section of the right diameter. The various sections of the connection can be removed by gripping firmly and unscrewing (Fig. 7). If simultaneous use of two drainage tubes of any diameter is necessary, join the single connections of the drainages to the suction

tube by means of the double Y-shaped connection.

3. Take care to shut the clamp.

4. With the bellows still compressed, remove the plug and connect up the suction tube connection by screwing it into the hole on the top of the bellows.

The clamp must be kept closed at this stage.

5. Hook the spring eyelet to the fastening by the hole on top, then to the fastening on the bottom.

6. Start suction by opening the clamp to maximum. Hang Emodren from the bed using the plastic suspender clip. When suction is complete, remove the drainage system.

Replacing the bellows.

When the system is in a fully extended position, it has completed its suction action. If it's necessary to continue, the bellows has to be replaced with a new one, as follows:

- close the clamp;
- replace the bellows using the dilatation spring;
- resume suction by opening the clamp.

The used bellows, containing the drained liquids, must be closed with its plug to allow safe disposal.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Re-use of single devices may involve the following risks:
 - deterioration of the performance of the devices;
 - infection or contamination by previous use on another patient.
- The device must only be used by qualified medical staff.
- After use, the device must be treated and disposed as a hospital waste, taking all the necessary precautions and as required by current provisions of law. Sharp parts must be disposed of in special rigid containers.
- The bellows is graduated to check the quantity of blood suctioned. The quantity in ml must be considered approximate.
- The use of original components is recommended. Medicalplastic declines all responsibility for harm caused to patients by the use of 3rd party products or use other than those specifically prescribed.
- The device is supplied sterile. Do not use after the expiry date or if the packaging has been opened or damaged.
- In case of a serious incident, please promptly notify the manufacturer and your Competent Authority.
- The device has not been specifically tested on pregnant and lactating women and paediatric patients. Any use on these populations should be evaluated by the physician.