emodren®

MD


Fig. 1



Fig. 2

Medicalplastic s.r.l.
Via Saverio Mercadante, 15 - 20124 Milano - Italia
tel. 026704424 - Email: info@medicalplastic.it - www.medicalplastic.it
MADE IN ITALY

500905

Rev.: 0 25/11/2024

IT MANTICI PER SISTEMA EMODREN

DESTINAZIONE D'USO

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EN BELLOWS FOR EMODREN SYSTEM

INTENDED USE

Bellows are a component of EMODREN, a continuous suction drainage system of the surgical wound to prevent post-operative haematomas.

DESCRIPTION

Bellows create suction and allow the collection of drained liquids.

Replacement bellows are available in the following models :

- 1) Bellows for Emodren 500ml Code 10012581 (**Fig. 1**) and
- 2) Bellows for Emodren Easy Code 40050030 (**Fig. 2**).

INSTRUCTIONS FOR USE

Figs. 1-2 Composition of bellows

- A. BELLOW**
- B. PLUG**
- C. SUSPENDER CLIP**

If it is necessary to continue the suction action of the system, the used bellows must be replaced with a spare one by proceeding as follows:

Bellows for Emodren 500 ml

- Close the clamp
- Release the bellows from the suction tube connection of the set in use.
- To reconnect the replacement bellows to the system, proceed as follows :
 - The bellows should be kept compressed, as in the sales package, and should be connected to the suction tube connection. Re-use the spring recovered from the system by first connecting it to the upper hole and then to the lower hole.
- Resume suction by opening the clamp.
- The used bellows, containing the drained liquids, must be closed with the plug to allow safe disposal.

Bellows for Emodren Easy

- Close the clamp
- Release the bellows from the suction tube connection of the set in use.
- To reconnect the replacement bellows to the system, proceed as follows :
 - Manually compress the bellows and, keeping it compressed, connect it to the suction tube connection by screwing it into the upper hole
- Resume suction by opening the clamp
- The used bellows, containing the drained liquids, must be closed with the plug to allow safe disposal.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

- Re-use of single devices may involve the following risks:
 - deterioration of the performance of the devices;
 - infection or contamination by previous use on another patient.
- The device must only be used by qualified medical staff.
- After use, the device must be treated and disposed as a hospital waste, taking all the necessary precautions and as required by current provisions of law.
- The bellows is graduated to check the quantity of blood suctioned. The quantity in ml must be considered approximate.
- The bellows are intended for use exclusively with drainage sets manufactured by Medicalplastic. We accept no liability for any harm caused to patients by the use of 3rd party products or use other than those specifically prescribed.
- The device is supplied sterile. Do not use after the expiry date or if the packaging has been opened or damaged.
- In case of a serious incident, please promptly notify the manufacturer and your Competent Authority.
- The device has not been specifically tested on pregnant and lactating women and paediatric patients. Any use on these populations should be evaluated by the physician.